

PRESCRIPTION

PANTHERA X3


Patient I	D:		
Dentist:			

Case identifier #:

TYPE OF BITE PROVIDED


■ MAXIMUM

The appliance will be set to -4mm of the maximum protrusion to allow incremental protrusion up to the maximum. Not recommended if maximum protrusion is less than 5mm.



STARTING POINT

The appliance will be set to start at the provided bite position, to allow incremental forward movement.



VERTICAL DIMENSION

Close or open to optimise the device

☐ Keep it, call if major changes needed

IS THERE A LATERAL DEVIATION DURING PROTRUSION?

Yes No

ELASTIC NOTCHES

Elastic notches delivered with all X3 devices.

FRAGILE TEETH, **CROWNS AND OR PONTICS:**

Tooth #:

BIOMATCH (use optimal values):

Yes No

* If YES checked, skip to section 5.

UPPER PLATEAU













UPPER BAND


















ANTERIOR WITH CONTACT (!)







LOWER BAND

















EXTRA OPTIONS Prefer distal wrap Do not cover 3RD molars Upper

Lower

☐ Add a Panthera morning repositioner (additional cost)

COMMENTS

S	IGN	TAI	'UR	RΕ

☐ Do not call me if design changes are needed.

X